# Aquamatic™ II

- D** Gebrauchsanweisung
- UK** Operating instructions
- F** Mode d'emploi

# Inhalt

Seite

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## Table of contents

Page

<b>1.0</b>	<b>Introduction</b>	
1.1	Notes on operating instructions.....	5
1.2	Intended use .....	7
1.3	Function .....	7
1.4	Explanation of symbols .....	9
<b>2.0</b>	<b>Safety advice</b>	<b>12</b>
2.1	Advice for the user for maintaining hygienic status.....	13
<b>3.0</b>	<b>Installation and start-up</b>	
3.1	Standard delivery .....	17
3.2	Illustrations.....	18 - 21
3.3	Connections .....	23
3.3.1	Electrical connection .....	23
3.3.2	Connecting a nystagmograph .....	23
3.3.3	Connecting a potential balance.....	23
3.3.4	Water connection .....	23
3.3.5	Connecting the irrigation handle .....	25
3.4	Start-up .....	25
<b>4.0</b>	<b>Operation</b>	
4.1	Adjusting temperatures .....	27
4.2	Selecting temperature levels.....	27
4.3	Adjusting stimulation time .....	27
4.4	Adjusting flow quantity .....	27
4.5	Description of operating modes .....	29
4.5.1	Thermal process for reducing bacterial count.....	29
4.5.2	Irrigation mode .....	29
4.5.3	Stimulation mode .....	31
4.5.4	Energy saving mode .....	31
<b>5.0</b>	<b>Cleaning</b>	
5.1	General information on cleaning and disinfection .....	33
5.2	Recommended instrument disinfectants .....	35
5.3	Recommended surface disinfectants .....	35
<b>6.0</b>	<b>Maintenance and service</b>	
6.1	Functional test.....	40
6.2	Visual inspection .....	41
6.3	Electrical test.....	41
6.4	Decalcification.....	41
6.5	Replacement of filter .....	42
<b>7.0</b>	<b>Trouble shooting</b>	<b>47</b>
<b>8.0</b>	<b>Accessories and spare parts</b>	
8.1	Accessories.....	50
8.2	Spare parts.....	50
<b>9.0</b>	<b>Specifications</b>	<b>53</b>
<b>10.0</b>	<b>Disposal</b>	<b>55</b>

General Standard Terms and Conditions

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## 1.0 Introduction

### 1.1 Notes on operating instructions

- The operating instructions contain important information on how to operate the Aquamatic™ II safely, correctly and effectively. The operating instructions are intended for new operating personnel and also should be continued to be used as a reference manual. Reading and understanding the operating instructions will assist the user in avoiding risks, reducing repair costs, lowering down-time, increasing reliability and extending the service-life of the equipment. **The operating instructions must always be kept available and near the equipment.** Prior to first use, please read the chapter "Safety advice", to be prepared for any possible dangerous situations.

The basic principles are:

**Judicious and careful work provides the best protection against accidents!**

Operational safety and readiness for use depend not only on your capabilities, but also on **care and maintenance** given to the Aquamatic™ II. Regular cleaning and service work are a must. Major maintenance and repair work should be carried out only by expert personnel authorised by producer. For repairs, only use the Producer spare parts to preserve the equipment warranty and ensure operational safety, readiness for work and to maintain the value of the equipment.

- The product Aquamatic™ II bears CE Marking CE-0124 according to the EEC guideline of the council for medical products 93/42/EEC and meets the basic requirements of this guideline.
- The quality management system applied at producer has been certified according to international standards EN ISO 9001 and EN 13485.
- Reprints, and text extracts are allowed only with written permission from Producer.

#### Abbreviations / symbols in these operating instructions:

- Indicating a list
  - Subdivision of a list/activity

The recommended sequence must be followed in each case!

Indicating particularly important advice!

↪ Describing the effect of an activity.

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UK

## 1.2 Intended use

- Compact unit for stimulation of the vestibule. Distinguished by its easy operation, it can be used in a hospital and in a practice.
- May also be used for irrigations of the auditory canal during examinations.
- The unit is available with different water flow quantities for the stimulation mode resp. irrigation mode. The respective flow quantity of the special unit is stated on the type plate.
- This unit may not be used if the auditory canal is inflamed or contaminated.
- This unit is designed for short-time use on patients.

## 1.3 Function

- Operating the main switch (❶, fig. 1) automatically activates the thermal process for reducing bacterial count (see section 4.4.1).
- The unit automatically switches to the irrigation mode where the auditory canal can be rinsed with water warmed up to 37°C.
- The unit operation can be changed to the stimulation mode for stimulating the vestibule with a reduced flow of water. The Aquamatic™ II is equipped with a timer for preselecting the stimulation time.

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






UK

## 1.4 Explanation of symbols

	Pay attention to operating instructions ! as to ISO /7000/0434 DIN 30600/1008 IEC 348
	Type B equipment as to IEC 417
	Fuse as to IEC 417/5016, DIN 30600/0186
°C	Temperature in degree centigrade
S	Timer adjustment in seconds
	Start
	Stop
	Timer
	Cold stimulation level
	Warm stimulation level
	Irrigation level (water at a temperature value of 37°C)
	Heating ON
	Heating OFF (energy saving mode)
	Control output for connecting a nystagmograph (graphical recorder as to DIN 30600, IEC 417 5192)
	Equipotentiality DIN 30600 495, ISO 417 5021
	Connection for waste water
	Connection for water
Low	Reduced flow for stimulation of the vestibule
High	High flow for rinsing the auditory canal

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⇒ Page 14 - 15

## 2.0 Safety advice

- The Aquamatic™ II unit is produced according to IEC 601/EN 60601 and listed in the following classes:
  - VDE Class of protection 1
  - Class IIa (EEC 93/42).
- The equipment must only be connected to a power outlet with a correctly installed ground/earthing contact according to local and state regulations.
- The unit must be installed by a specialist authorized by GN Otometrics (see chapter 3.3).
- The Aquamatic™ II may be used in **supervised operation** by qualified personnel which have been authorised by GN Otometrics and which have been trained in the operation of the equipment (IEC 601-1/EN 60601-1).
- The mains voltage specified on the type plate must match the power supply system.
- Prior to every use of the equipment, verify that the device is technically safe and in proper working condition. **Damaged leads** must be replaced immediately!
- Correct configuration in assembly of country-specific connections
 

green/yellow:	protective conductor (PE)
blue:	neutral conductor (N)
black resp. brown:	phase (L)
- The control panel must be visible and in reach of the operator. Ensure the equipment is stable prior to use.
- In the thermal process for reducing bacterial count hot water is being ducted through the irrigation handle. Please do **not** take the handle from its holder or spray water!
- Metal parts can be hot!
- Prior to spraying, the water temperature must be checked by the user!
- Switch off main power switch after use and **close the water valve** of water supply.
- The Aquamatic™ II should only be operated in rooms used for medical purposes, and never in areas subject to explosion hazards.
- All additional equipment, which is connected to the analog and digital interfaces of the unit, must meet the requirements of relevant EN specifications (for inst. EN 60950 for data processing equipment and EN 60601 for electrical medical appliances). In addition, configurations must satisfy system specification EN 60601-1-1. When additional equipment is connected to the signal input or signal output section on the unit, the person making the connection is deemed "a system configuration operator" and as such is responsible for meeting the requirements of system specification EN 60601-1-1. For answers to additional questions, please contact your local specialist supplier or GN Otometrics technical service.
- The irrigation jet must not get into contact with contaminated material.
- Only use the warm water jet when a tube tip is installed!
- Take caution to avoid injury to the eardrum when introducing the water jet !
- For hygienic reasons jet connections must be changed after each patient. This also prevents retrograde germination of the warm water system.
- Use only for irrigation of auditory canal !
- GN Otometrics is not liable for personal injury and damage to property if
  - original GN Otometrics parts are not used,
  - the operating instructions are not completely followed,
  - assembly, new settings, alterations, extensions and repairs have been carried out by personnel not authorised by GN Otometrics.

## 2.1 Important instructions for the maintenance of the hygiene status of warm water irrigation units

For loosening cerumen in the auditory meatus and for the stimulation of the labyrinthine system, GN Otometrics offers the warm water units Hygrotherm plus (37°C) and Aquamatic™ II (20° - 47°C).

These units heat the drinking water which comes from the household connection to the preselected temperature. The operating instructions must be followed, or the increased number of pathogens in the rinse water can affect the health of sensitive patients.

### Requirement for taking into use

- ☐ The water which is provided by the household connection must at least meet the WHO guidelines or the country-specific guidelines for drinking water.

### Connection

- ☐ Fresh water must flow through the wall connection for a period of approx. 1 minute prior to the unit being connected.\*
- ☐ Ensure that the installation is hygienically carried out! Disinfect the threaded connections with 70 % alcohol\* before screwing them into place.
- ☐ Wear disposable gloves\* when carrying out maintenance work on parts which come into contact with the water.
- ☐ Carry out thermal cleaning at least 3 times before patient use.\*

\* This work is carried out by the GN Otometrics service technician.

### Running operation

- ☐ Before beginning the surgical operation, switch on the unit and await the thermal cleaning process.
- ☐ After the device had been idle for a long period of time, (weekend, during holidays, etc.) the thermal cleaning process must be initiated and repeated.
- ☐ We recommend that the device be switched off and on every 2 hours and not longer than 4 hours to start the thermal cleaning mode.
- ☐ Immediately replace and autoclave contact parts such as the jet connection each time it is used for a different patient (for prevention of retrograde germification).
- ☐ Before every use, spray water and check the temperature.

### Service

- ☐ Ensure that the suspended particle filter is hygienically changed (see commissioning) (dispose of filter and carefully clean the filter glass). Ensure that no contact is made with the contaminated parts.
- ☐ Carry out thermal cleaning after each service.
- ☐ We recommend that the number of germs in the water be regularly determined at the input and output of the unit.

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### 3.1 Standard delivery

#### Quantity

1	Aquamatic™ II basic unit
1	Irrigation handle
1	Double-barrelled hose
3	Jet connection short, straight (80 mm)
1	Splash protection
1	Rubber bush for splash protection
2	Three-holes gasket
1	Hose tips for jet connection (30 Stk)
1	G3/4i-G1/4a-screwing (water-tap G3/4 on filter G1/4)
1	Cap nut G3/4
1	Water filter complete
2	Gasket (f.G1/4a) 13x18x1
1	G1/4a-G3/4a-reducer (filter G1/4 on water hose G3/4)
1	Water supply hose G3/4i, L = 3 m
1	Waste water hose G1/4i, L = 3 m
2	Gasket (f. G1/4i) 10x15x1
1	Cardboard 455 x 340 x 250 mm
1	Padding for cardboard
1	Operating instructions

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D

Fig. 1. Front view


- ① Main switch  
② Key switch f. heating ON/OFF (energy saving mode)  
③ Coding switch for warm-stimulation level  
④ Coding switch for cold-stimulation level  
⑤ Coding switch for stimulation time  
⑥ Key switch for selecting the irrigation level (37°C)  
⑦ Key switch for selecting the warm-stimulation level (e.g. 44°C)  
⑧ Key switch for selecting the cold-stimulation level (e.g. 30°C)  
⑨ Key switch for start/stop of the stimulation  
⑩ Temperature display (two-figure number, resolution 1°C), current value indication  
⑪ Display of stimulation time (two-figure number, resolution 1s)

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**D****Fig. 2. Rear view**

- ❶ Type plate
- ❷ Control output for a nystagmograph
- ❸ Unit plug with fuse compartment
- ❹ Equipotentiality connection

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D

Fig. 3. Side view


- ① Support for irrigation handle
- ② Connection for double-barrelled hose
- ③ Connection for waste water hose
- ④ Connection for water supply hose

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D

Fig. 4. Irrigation handle

- ① Handle
- ② Connection for double-barrelled hose
- ③ Three-holes gasket
- ④ Knurled screw
- ⑤ Release lever
- ⑥ Valve plunger
- ⑦ Jet connection
- ⑧ Splash protection
- ⑨ Cap
- ⑩ Double-barrelled hose
- ⑪ Lever for water flow

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### 3.3 Connections

#### 3.3.1 Electrical connection

- Connect the power cable to the socket on the unit (❸, fig. 2).
- Insert the power plug in a correctly installed socket with a ground/earthing contact.

#### 3.3.2 Connecting a nystagmograph

- When controlling an ENG (electro-nystagmograph) or a CNG (computer-nystagmograph) at output (❷, fig. 2), please connect only recording equipment approved by GN Otometrics. Connecting cable available from A GN Otometrics (see section 8.0).
- At the end of the stimulation time, a trigger signal for a nystagmograph is available at the rear 3-pole DIN socket. This outlet is an electrically isolated photo transistor. The emitter connection which has to be connected to the ground line of the trigger inlet at the nystagmograph is attached to pin 3. The open-collector connection is located on pin 1 and is to be connected to the positive supply voltage (+5 V or +12 V) by means of a pull-up resistance. The maximum collector voltage of the transistor may not exceed 80 MA.

#### 3.3.3 Equipotentiality connection

- Connection for potential compensation (❹, fig. 2). Connecting cable available from GN Otometrics (see section 8.0).

#### 3.3.4 Water connection

- Local requirements:
  - Water-tap with G3/4" external thread.
  - Drinking-water!
  - Input water pressure: +2000...+6000 hPa.
  - Feed temperature: +7...+19.5°C (recommended), but at least 0.5°C below the required lowest cold stimulation value.
  - Drainage fittings (Art.No. 502.0880.0) as to DIN 1988.
- Connection to water supply and waste water:
  - Prior to connecting the unit to the water supply the feed line must be flushed clean by opening the water tap for a minute and allowing free flow of the tap water.
  - Connect the filter unit to the water-tap by means of adaptor and gasket.
  - Insert the gaskets in the cap nuts of the water supply hose.
  - Join the nuts with the filter connection and the unit connection (❹, fig. 3).
  - Insert the gaskets in the cap nuts of the waste water hose
  - Join the nuts with the unit connection (❸, fig. 3) and the drainage fittings.

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UK

### 3.3.5 Connecting the irrigation handle

- Insert a "three-holes gasket" (③, fig. 4) in the double-barrelled hose and screw it (⑩, fig. 4) on the unit connection (②, fig. 3).

Use only "three-holes gaskets" or the unit will not function correctly!

- Insert a "three-holes gasket" in the irrigation handle and screw it on the free end of the double-barrelled hose.
- Open the water-tap.  
Check whether all connections are tight.

### 3.4 Start-up

- Insert the handle in its holder; the jet must point to the reverse side of the unit.
- Switch on the unit (①, fig. 1).
- The unit will display a test mode with digital numbers "8 8" and an audible signal.
- Automatic activation of the thermal process will begin for reducing bacterial count.  
Duration: 5 min.
- When the thermal process is complete, the unit will switch to the irrigation mode (37°C).
- When the irrigation mode is not used (handle is not removed from its support) the unit change into the energy saving mode after 5 min.

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
- ① Stellhebel
- ② Justierschraube

Fig. 6. Irrigation handle

- ① Lever for water flow
- ② Adjusting screw

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UK

## 4.0 Operation

- Before switching on the unit, make sure that the handle is in its holder (jet must point to reverse side of the unit).

### 4.1 Adjusting temperatures

- Number of temperature levels: 3
  - One level fixed to 37°C irrigation temperature.
  - Two variable temperature levels (20°C - 47°C) (48°C + 49°C only for testing purposes).
- Temperate setting by coding switches (③, ④, fig. 1).
  - left switch: for adjusting the "ten" partition
  - right switch: for adjusting the "one" partition
    - ↳ lower keys (+): Temperature increase
    - ↳ upper keys (-): Temperature decrease
- Standard settings:
  - Level for irrigation mode: 37°C fixed
  - Level for cold stimulation: 30°C
  - Level for warm stimulation: 44°C

### 4.2 Selecting temperature levels

- After operating the main power switch, the automatic activation of the thermal process begins for reducing bacterial count, followed by the automatic activation of the temperature level "irrigation mode" (37°C).
- For selecting the desired temperature level, use the respective key (⑥, ⑦, ⑧, fig. 1).
  - ↳ Display of the active level by illuminated diodes.
  - ↳ Display of the water temperature (current value) in °C.

### 4.3 Adjusting stimulation time

- By means of coding switch (⑤, fig. 1).

### 4.4 Adjusting flow quantity

- With the lever for water flow (①, fig. 6) you can adjust the flow quantity at the irrigation handle.
- 2 adjusting possibilities:
  - High:** Flow quantity for rinsing the auditory canal
  - Low:** Flow quantity for stimulating the vestibular organ
- With the adjusting screw (②, fig. 6) you can precisely pre-select the flow quantity.

Periodically check the water flow for the stimulation of the vestibular organ and readjust, if necessary.

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UK

## 4.5 Description of modes

### 4.5.1 Thermal process for reducing bacterial count

#### Purpose:

Reduction of bacteria count in the warm water circuit to prevent or clear a contaminated system.

#### Activation:

Occurs each time the unit is switched on.

- Condition for starting the correct thermal process for reducing bacterial count:
  - The irrigation handle must be inserted in the holder on the side (jet pointing to reverse side of the unit).
    - ↳ If not, a warning signal is given (warning also sounds when the handle is removed during the thermal cleaning procedure).
    - ↳ Additionally, the error code "F0" is shown (disappears when handle is returned).
  - In the thermal process for reducing bacterial count, hot water is passing through the irrigation handle. Please do **not** take handle from its holder or splash off water!
  - Metal parts of the irrigation handle and of the double-barrelled hose will be hot! Please do not touch metal parts during the thermal cleaning procedure!
  - During the thermal process for reducing bacterial count, the temperature display shows a 'd' on the first segment and the remaining time (in minutes) on the second segment.
- After cooling down to 37°C the unit will automatically switch into irrigation mode.

### 4.5.2 Irrigation mode

#### Purpose:

Rinsing of the auditory canal.

#### Properties:

- Temperature: 37°C, permanently set.
- Position of the adjusting lever: **High** (see fig. 6, page 26).

#### Activation:

- Automatically after completion of thermal cleaning procedure and taking out the handle or by key (ⓐ, fig. 1).
- When the irrigation mode is not used (handle remains in holder) the Aquamatic™ II changes into the energy saving mode after 5 min.
- Reactivation of the irrigation mode by taking the handle out of the holder or activating a key on the temperature control unit.
- Repeated activation of the 37°C key effects switching-off of the heating system.
  - ↳ Water with feed temperature is available.

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Stimulation of the vestibule.

#### Properties:

- Temperature: according to pre-set cold or warm water stimulation level
- Position of the adjusting lever: **Low** (see fig. 6, page 26)
- Duration: as pre-set by means of timer.

#### Activation:

- Initial selection of the type of stimulation is done by actuating either warm-stimulation or cold-stimulation key (7, 8, fig. 1) (see section 4.1 for pre-setting of temperature).  
↳ Water with pre-set temperature value flows to handle.

- Wait for 20 seconds.

- Operate the "timer-start-key".

- Preparation for stimulation:
  - Operate unit 15 seconds in holding position, so that the jet can be positioned correctly in the auditory canal.

During this period, the water is directed, via the handle, into the drain in order to guarantee an optimum temperature stability; with an optical indication by the flashing figure in temperature display.

Do not start stimulation by operating the release lever at the handle until you hear the audible signal!

- Carrying out the thermo-stimulation by actuating release lever at handle (5, fig. 4).
- At the end of the stimulation period a control signal for a recording unit is issued at the nystagmograph output.
- After completion of the stimulation the unit remains in a holding position for 15 seconds so that the handle may be removed from the auditory canal.  
↳ During this time, water flows into the drain; again an optical indication is made by flashing temperature value.
- Repeated actuation of the "timer-start-key" leads to stopping the timer operation.
- Second actuation of the currently active key effects deactivation of the pertaining stage.  
↳ Heating is switched off completely.  
↳ Stimulation with cold water (feed temperature).

### 4.5.4 Energy saving mode

#### Purpose:

- Reduction of energy consumption.

#### Activation:

- Actuation of key "heating on/off" (2, fig. 1).  
↳ Heating is switched off.
- After five minutes of non-use of the unit.

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UK

## 5.0 Cleaning

### 5.1 General information on cleaning and disinfection

After use, all parts which come into contact with the patient (rinsing connection, jet connection, hose tips) must be removed and disinfected resp. autoclaved! You may use all cleansers and disinfectants listed in section 5.2.

The surfaces of the Aquamatic™ II resist most of the common surface disinfectants.

#### However, do not use

- disinfectants containing concentrated organic or anorganic acids or bases, since these may cause corrosion damage;
- disinfectants containing chloramides, phenol derivatives or anionic tensides, since these may cause stress cracks in the plastic material used for the housing of the unit.

For cleaning and disinfection, you may also use disinfectant sprays or disinfectant tissues.

Switch off the unit by the main power switch before starting cleaning and disinfection.

Wipe the unit surface with a cloth moistened with a cleaning or disinfecting solution. Make sure that no liquid penetrates the unit. You may use all cleansers and disinfectants listed in section 5.3.

Always observe the instructions for use by the manufacturer of the disinfectants, above all, concentration specifications.

The above measures for cleaning and sterilizing do not replace the respective precautions for operating the unit.

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F

⇒ Page 36

## 5.2 Recommended instrument disinfectants

Disinfectant	Ingredients	(in 100 g)	Manufacturer
GIGASEPT FF (Application concentrate)	succinic acid dialdehyde dimethoxy tetrahydrofurane corrosion protection components non-ionic tensides and odoriphores	11.0 g 3.0 g	Schülke & Mayr, Norderstedt
Sekusept PLUS (Application concentrate)	glucoprotamine non-ionical tensides dissolvents, complexing agents	25.0 g	Henkel, Düsseldorf
Mucozit-T new (Application concentrate)	bi3 (3-aminoprophyl)laurylamine alkyldimethylbenzylammoniumchloride cocospropyldiamine-1,5 guanidinium-acet.	8.0 % 19.0 % 7.0 %	Merz & Co., Frankfurt/Main

## 5.3 Recommended surface disinfectants

Disinfectant	Ingredients	(in 100 g)	Manufacturer
TERRALIN (Application concentrate)	Benzal conium chloride Phenoxypropanole	20 g 35 g	Schülke & Mayr, Norderstedt
QUATOHEX (Application concentrate)	Didecyl dimethyl ammonium chloride Benzal conium chloride Bi-guanidinium acetate Polymer biguanid Cleaning agents	14 g 10 g 7,5 g 0,5 g	Braun, Melsungen
Incidin Plus (Application concentrate)	Glucoprotamin Nonionic tensides Solvents, complexing agents	26,0 g	Henkel, Düsseldorf
Pursept-A (Disinfectant spray or disinfectant cloths)	Ethanol Glyoxal QAV	38,9 g 0,1 g 0,05 g	Merz & Co., Frankfurt/M.

If using aldehyde-containing or amine-containing disinfectants at the same object, this may result in discolourations.

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⇒ Page 43 - 45

## 6.0 Maintenance and service

Decalcification and changing of filter (see sections 6.4, 6.5) should be carried out immediately by the user when required.

### Safety-technical inspection

#### 6.1 Functional test

##### Test

- ☐ ● Switch unit several times on and off.
  - ↳ Correct function of mains pilot lamp and solenoid valve (clacking).
  - ↳ Automatic activation of the thermal process for reducing bacterial count.
- ☐ ● Shut-off main water tap and observe temperature indication.
  - ↳ After a few seconds error signal will appear "F1" (= no water).
- ☐ ● Testing the displays.
  - ↳ All optical and audible indicators issue their signal during the first second after switching on (in particular segmental test "8 8 ").
- ☐ ● Set all temperature setting switches on 40°C and test consecutively by touching the corresponding keys whether in all ranges 40°C is reached in ACTUAL temperature indication.
  - Spray off and check water temperature by hand.
  - Carry out the same test with temperatures 37°C and 28°C.
  - ↳ The particular combination of temperatures 40°C, 37°C and 28°C, as well as the test in all ranges ensures a comprehensive test of all contacts of the pre-setting switches, in ON- as well as in OFF-position, and also testing of all components involved in temperature control of the DESIRED temperature control system. Therefore, carrying out a test in one range only, or at only one temperature is not permissible.
- ☐ ● Set temperature at 49°C.
  - ↳ The electronic excess-temperature protection must respond (in the temperature display error code "F7" appears).
- ☐ ● Measure temperature of water splashed off at standard settings: 30°, 37° and 44°C.
  - ↳ The value may not deviate more than  $\pm 0,5^\circ \pm 1/2$  digit, i.e. 1°C in total. Typically, however, values are better than  $\pm 0,5^\circ\text{C}$ .

##### **This measurement is critical !**

Measuring probes or thermometers which are too large may cause faulty measurements because of substantial heat removal, inertia, and because their temperature sensitivity is not localised on a point but spread over a larger surface. Use measuring instruments with an immersion probe, having a minimum precision of 0.3°C and a resolution of 0.1°C.

After setting the temperature, wait at least 1 minute before measuring to allow temperature stabilization of the double-barrelled hose and handle.

See also section 4.5.3!

## 6.2 Visual inspection

### Test

- ☐ ● Check unit for external damage (damage to cables, hoses, filter and irrigation handle).
- ☐ ● Check hose system for external indications of leaks (drops, water spots, corrosion deposits etc.).
- ☐ ● Inspect unit, disconnected from the main power supply and with open housing, for internal damage, leaks, corrosion, and deposits (mineral scale). Close housing again and tighten screws. Re-connect main power cable.
- ☐ ● Check fuses for correct protection value and melting characteristics.
- ☐ ● Check sign plates and inscriptions for completeness and readability.


## 6.3 Electrical test

### Test

- ☐ ● Measure **protective conductor resistance** according to VDE 0751 (max. 0.3 Ohm).
- ☐ ● Measure and compare with "first value measured" of **unit leakage current** according to VDE 0751 (measured during operation).
  - ↳ The "value measured first" may not be exceeded by more than 50%.
- ☐ ● Document and record the evaluation and tests according to VDE 0751.

## 6.4 Decalcification

- ☐ ● If the tap water in the area, where the unit is used, is known for mineral precipitation, preventive decalcification should be carried out:
  - Close water tap.
  - If applicable, dismantle filter when mounted vertically at water tap.
  - Unscrew screwed glass of filter.
  - Fill screwed glass with approx. 100 ml of decalcifier on acetic / citric, or lactic acid base (for inst. Citrosteril, by Fresenius).
  - Screw glass on again.
  - Install filter vertically again, if applicable.
  - Open water tap again.
  - Switch on the unit.
  - ↳ Decalcification takes place automatically in the thermal process for reducing bacterial count. After the unit has been in this mode for 5 minutes the water passages have been decalcified and cleaned, and also sufficiently flushed. The Aquamatic <sup>TM</sup> II is again fully ready for operation. When water passages are heavily calcified, a modified procedure may become necessary, with a prolonged reaction period for the decalcifier.
- Recommendation: Calcification safety device 502.0995.0.

*Fig. 7. Water filter*

- ❶ Filter body
- ❷ Filter cartridge
- ❸ Sealing ring
- ❹ Plexiglass cup

## 6.5 Replacement of filter

- When the filter cartridge is very dirty (dark colour) the filter cartridge is to be exchanged as follows:
  - Close water tap of supply.
  - Switch on unit to reduce pressure in the filter body.
  - Switch off unit after appearance of "F 1".
  - Unscrew Plexiglass cup (❹, fig. 7) from filter body (❶, fig. 7) by turning anti-clockwise.
  - Extract filter cartridge (❷, fig. 7) and replace by new one (Art.No. 502.0891.0).
- The filter element must not come into contact with contaminated objects in order to prevent the ingress of germs into the water system.

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


## 7.0 Troubleshooting

If, errors cannot be corrected inform the GN Otometrics service staff of your supplier. Do not attempt to repair the unit

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## 8.0 Accessories and spare parts

### 8.1 Accessories

<u>Description</u>	<b>REF</b>
Calcification safety device.....	502.0995.0
Jet connection extra long, straight (110 mm).....	508.0429.0
Rinsing connection (Teflon) for stimulating the Vestibular organ .....	501.0355.0
Extension of supply hose G3/4a-G3/4i, L =1,5 m .....	501.0315.6
Waste water hose, special length, G1/4a-G1/4a (please indicate desired length) .....	502.0881.1
Drainage fittings as to DIN 1988 .....	502.0880.0
Connecting cable for a nystagmograph .....	502.0850.0
Connecting cable for the equipotential compensation, L = 5m .....	008.0596.0

### 8.2 Spare parts

<u>Description</u>	<b>REF</b>
Aquamatic™ II, basic unit, 220-240 V AC, 50/60 Hz.....	
Aquamatic™ II, basic unit, 110-127 V AC, 50/60 Hz .....	
Irrigation handle 4+ .....	502.0963.0
Double-barrelled hose.....	510.0412.0
Jet connection short, straight (80mm).....	508.0427.0
Splash protection .....	501.0331.0
Rubber bush for splash protection.....	501.0331.1
Throttling nozzle 0.5 mm (150 ml/min).....	502.0946.0
Throttling nozzle 0.7 mm (250 ml/min).....	502.0946.1
Three-holes gasket .....	501.0320.0
Hose tips for jet connection (30 pcs.).....	502.0844.0
G3/4a-G3/4i, L = 3 m .....	502.0768.1
Cap nut G3/4 (to be included in order).....	052.0035.1
Water filter complete .....	502.0890.0
Filter element for water filter .....	502.0891.0
Gasket (for G1/4a) 13x18x1.....	055.0018.0
G1/4a-G3/4a reducer (filter G1/4 to water hose G3/4).....	502.0724.1
Supply hose with connection G3/4i, L = 3 m.....	501.0315.7
Waste water hose with connection G1/4, L = 3 m.....	502.0882.0
Gasket (für G1/4i) 10x15x1 .....	055.0003.0
Power cable with protective contact socket (Germany), L = 3m .....	507.0859.0
Power cable with British plug BS1363A (UK), L = 2,5m.....	008.0762.0
Fuse F 8 A (quick-acting) 250 V breaking capacity H (1500 A).....	008.0767.0
Fuse for special voltage (110-127VAC): F 15A (quick-acting) breaking capacity 10kA .....	008.0766.0
Front foil Aquamatic™ II .....	
Operating instructions.....	502.0900.i

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F

⇒ Page 54

Voltage range	220 - 240 V~ ± 10 % 50/60 Hz Special voltage: 100 - 127 V~ ± 10 % 50/60 Hz
Current input	6.8 A (220 - 240 V~); 14.8 A (100 - 127 V~)
Power cable	3 m
Power consumption	1500 W
Operating modes	Thermal process for reducing bacterial count; Irrigation mode; Stimulation mode; Heating off (energy saving)
Stimulation time	adjustable by timer from 1 up to 99 sec.
Timer indication	Accuracy ± 0.5 s ± ½ digit
Water temperature	20°C - 47°C
Temperature indication	Accuracy ± 0.6 °C ± ½ digit (can only be guarantee with unchanged ambient conditions)
Feed temperature of the water	+7...+19.5°C
Water quality	Drinking-water
Water flow in the irrigation mode (high flow)	450 ml/min: standard unit and variant 250/450 ml/min 400 ml/min: variant 400/400 ml/min 500 ml/min: variant 500/500 ml/min
Water flow in the stimulation mode (low flow)	150 ml/min: standard unit 150/450 ml/min (adjustable) 250 ml/min: variant 250/450 ml/min (adjustable) 400 ml/min: variant 400/400 ml/min (adjustable) 500 ml/min: variant 500/500 ml/min (adjustable)
Operating pressure	at least: +2000 hPa (2 bar); max: +6000 hPa (6 bar)
Connections	Power connection ; control output for a nystagmograph; equipotential equalization; water supply G 3/4" external thread; water drainage G 1/4" external thread; connection for double-barrelled hose
Operating time	Continuous operation
Fuses	quick-acting 8 A (220-240 V~); quick-acting 15 A (110-127 V~)
Protective earth conductor resistance	0,1 Ω max.
Earth leakage current	0.5 mA max.
Enclosure leakage current	0.1 mA max.
Patient leakage current	0.1 mA max.
Ambient conditions	
Transport/storage	-20...+50°C; 5...90 % humidity without condensing, air pressure 700...1060 hPa
Operation	+10...+35°C; 20...80 % humidity without condensing, air pressure 700...1060 hPa
Dimensions HxWxD	14.5 x 37 x 31.5 cm
Weight	approx. 5.5 kg
Protection class (EN 60601-1)	I
Degree of protection	Type B 
Protection category	IPX0
Classification acc. to Annex IX EEC directions 93/42/EEC	Ila
CE marking	CE 0124
Rules applied	EN 60601-1: 1990 + A1:1993 + A2:1995 EN 60601-1-2: 1993 (EMV / EMC)
Canadian Classification	
Device group	Ear, nose, throat
PNC	77ETP
Risk Class	2
Description	STIMULATOR, CALORIC (WATER)

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UK

## 10.0 Disposal

- Packaging material, cardboard and/or PE foam, can be fully recycled or returned to your supplier.
- The Aquamatic™ II does not contain any hazardous materials.
- The material of the housing can be completely recycled.
- The component parts of the Aquamatic™ II must be correctly discarded and the materials should be carefully separated for recycling.
- Electronics circuit boards must be discarded using the appropriate recycling process.
- Used hose tips, which no longer can be disinfected, must be discarded into domestic waste immediately.

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This document is subject to change without notice.

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**EC - DECLARATION OF CONFORMITY  
FOR MEDICAL PRODUCTS**

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GmbH & Co. KG**